

Laser-aided plasma diagnostics of a helicon plasma source for fusion neutral beams

Ivo Furno

R. Agnello, Y. Andrebe, H. Arnichaud, P. Blanchard, A.A. Howling, R. Jacquier, G. Plyushchev

Swiss Plasma Center, Ecole Polytechnique Fédérale de Lausanne, Switzerland

Ph. Guittienne

Helyssen, Route de la Louche 31, CH-1092 Belmont-sur-Lausanne, Switzerland

M. Barbisan, R. Pasqualotto



Consorzio RFX, Corso Stati Uniti 4, I-35127 Padova, Italy

S. Bèchu

Université Grenoble-Alpes, CNRS, Grenoble INP, LPSC-INP23, 38000 Grenoble, France

I. Morgal and A. Simonin

CEA, IRFM, F-13108 St-Paul-lez-Durance, France




Access to H-mode regime

Achievement of a steady burn

Soft termination of discharge

ITER H&CD NBI	D-
Neutrals Power	16.5MW
Neutrals Energy	1MeV
Beam Current	40A
Pulse length	3600s
Neutralization eff.	~55%





H&CD NBI ¹	D-
Neutrals Power	16.8MW
Neutrals Energy	800KeV
Beam Current	34A
Pulse length	7200s
Neutralization eff.	>70%

[1] P. Sonato et al., Nucl. Fusion **57** (2017)

[2] A. Simonin et al., Nucl. Fusion **55** (2015)

Can helicon plasmas be used to produce negative ions for NBs?






6 independent
B field coils

$B_z < 800$ Gauss

water-cooled vacuum
vessel
(length 2m, diam.
0.4m)

Birdcage
resonant
antenna

Full diagnostic access




9 water-cooled copper legs

$D_{\text{int}} = 13 \text{ cm}$
 $L = 15 \text{ cm}$

Water OUT

Water IN


3 nF Mica Capacitors




Birdcage antenna

Optical emission spectroscopy


high throughput spectrometer





water-cooled end plate




OES and collisional-radiative code reveal negative ions in H₂/D₂ plasmas



- Principle of laser photo-detachment (PD)
- Cavity Ring Down Spectroscopy (CRDS)
- Langmuir probe (LP)-assisted laser photo-detachment
- Results
- Conclusions



S.J. Smith and D. S. Burch, Phys. Rev. Lett. 2, 165 (1959)



Principles of cavity ring down spectroscopy (CRDS)



$$I(t) = I_0 e^{-t/\tau_0}$$

$$t = \frac{L}{c(1 - R + ds n_{H^-})}$$


Neg. ion density

Negative ion rich shell

$$t_0 = \frac{L}{c(1 - R)}$$

Mirror reflectivity

$$n_{H^-} = \frac{1}{s} \frac{L}{cd} \left(\frac{1}{t} - \frac{1}{t_0} \right)$$



FWHM of laser pulse ~ 5ns

1

Avoid overlapping modes in the cavity

$$L_{\min} = c \Delta t / 2 = 75 \text{ cm} \quad L = 0.91 \text{ m}$$

Curv. radius of HR mirrors



2



Stability of the cavity


$$0 < \left(1 - L/r\right)^2 < 1 \quad r = 1 \text{ m}$$


Layertech HR mirrors with nominal $R = 99.994\%$ → $\tau_0 \sim 50 \mu\text{s}$

^[1] E. A. Bahaa, et al., Fundamentals of Photonics, John Wiley and Sons, Inc. (1991)










LP in e⁻ sat mode

$$i = i_{DC}$$




Negative ion rich shell

LP in e^- sat mode

$$i = i_{DC}$$



$$i = i_{DC} + i_{PD}$$




Negative ion rich shell


$$\frac{i_{DC}}{i_{PD}} = \frac{n_e}{n_{H^-}}$$

M. Bacal et al., RSI 2000

Diaphragm



IR mirrors



Beam dump


IR photod. for trigger

RF compensated L-shaped LP


1 Probe bias must be above the plasma potential


2 PD current independent from beam spot size

3 All neg. ions in the target volume must be PD






$$\frac{Dn_{H^-}}{n_{H^-}} = 1 - \exp\left(-\frac{sE}{hnS}\right)$$








R. Agnello et al., submitted to Nucl. Fusion





R. Agnello et al., submitted to Nucl. Fusion





- In RAID, helicon plasmas generate a significant negative ion population in H_2/D_2
- Successful use of CRDS and LP-assisted laser photo-detachment
- More laser-based diagnostics
 - ✓ Thomson scattering (P. Blanchard P57)
 - ✓ LIF
- Next: negative ion extraction

